

## W2100-System Battery Charging Set

**Charging Set for Storage Batteries for the Headstages of the W2100-System**  
**W2100-C-LC for with Batteries with Low Capacity (30 mAh)**  
**W2100-C-HC for Batteries with High Capacity (100, 200, 300 mAh)**

Two types of battery charger are available:

1. **W2100-C-HC** for storage batteries with high capacity (100, 200, 300 mAh)
2. **W2100-C-LC** for storage batteries with low capacity (30 mAh).

Do not mistake the charging devices! Do not recharge batteries without supervision!




### Applications

The charging device for the storage battery of the W2100-System headstage is USB powered.

To recharge the battery, please connect the charging set to an USB port of the computer and plug the battery via adapter cable into the connector of the charger.

The red LED on top indicates that the device is on. The yellow LED indicates that the storage battery is being recharged.

After approximately one hour the green LED indicates the end of the recharging process.



Dimension of the W2100 Charging Device:

58 mm x 20 mm x 10 mm

Weight of the charging device: 7 g

### Scope of delivery:

1 W2100-C-LC or W2100-C-HC charger and 1 adapter cable (W2100-C-B-Ch)

### W2100-C-HC:

W2100-B-100mAh-CA

W2100-B-200mAh-CA

W2100-B-300mAh-CA

### W2100-C-LC:

W2100-B-30mAh-CA